




**amplitude**




**eclipse**  
profil











Amorphous PLDLLA

Range Ø 7 to 10 mm

Progressive insertion




**Range :**

| Ø \ lgth | 7mm   | 8mm   | 9mm   | 10mm  |
|----------|---|---|---|---|
| 20mm     |   |   |   |   |
| 25mm     |  |  |  |  |
| 30mm     |   |  |  |  |

**Features :**

Progressive insertion  
5 branch star imprint



**Material :**


Polylactic acid PLDLLA



Screwing and tapping on Ratcheting handle



One piece : screwing and tapping



Screwing on 1.1 mm pin guide

- Stiff
- or
- Nitinol

| Designations  | References |
|---|------------|
| Polylactic interference screw ECLIPSE Profil Ø 7 lg 20mm  | 1-0402771  |
| Polylactic interference screw ECLIPSE Profil Ø 7 lg 25mm  | 1-0402772  |
| Polylactic interference screw ECLIPSE Profil Ø 8 lg 20mm  | 1-0402781  |
| Polylactic interference screw ECLIPSE Profil Ø 8 lg 25mm  | 1-0402782  |
| Polylactic interference screw ECLIPSE Profil Ø 8 lg 30mm  | 1-0402783  |
| Polylactic interference screw ECLIPSE Profil Ø 9 lg 20mm  | 1-0402791  |
| Polylactic interference screw ECLIPSE Profil Ø 9 lg 25mm  | 1-0402792  |
| Polylactic interference screw ECLIPSE Profil Ø 9 lg 30mm  | 1-0402793  |
| Polylactic interference screw ECLIPSE Profil Ø 10 lg 25mm | 1-0402702  |
| Polylactic interference screw ECLIPSE Profil Ø 10 lg 30mm | 1-0402703  |

Legal manufacturer : Biomatlante • 5 rue Édouard Belin • Z.A. Les Quatre Nations • 44360 • Vigneux de Bretagne • France

#### Service Clients – France :

Porte du Grand Lyon,  
01700 Neyron – France  
Tél. : **+33 (0)4 37 85 19 19**  
Fax : +33 (0)4 37 85 19 18  
E-mail : [amplitude@amplitude-ortho.com](mailto:amplitude@amplitude-ortho.com)

#### Customer Service – Export :

11, cours Jacques  
Offenbach, ZA Mozart 2,  
26000 Valence – France  
Tél. : **+33 (0)4 75 41 87 41**  
Fax : +33 (0)4 75 41 87 42