

Technical Sheet


SATURNE® II

Dual Mobility Cup
Cementless Ti / HA coated



 AMPLITUDE®

- Cylindro-spherical cup with 3mm cylindrical overhang
- Flat pole
- Equatorial press-fit: 1mm at diameter on all sizes
- Dual coating: titanium spray (80µm) and HA (80µm)



- Parallel concentric grooves on equatorial area
- Inverted grooves on same area



- Mirror-polished inner surface
- Rounded edges
- Material: Stainless steel M30NW

- Instrumentation available in straight or curved version
- One instrumentation tray



- Retentive liner for femoral heads Ø22.2 mm and Ø28 mm
- Large chamfer au at retentive threshold: later contact between cup and stem's neck, higher range of motion between stem and liner.
- Material: UHMWPE



References

CE 0086

DUAL MOBILITY CUP SATURNE® II cementless Ti / HA coated

Dual Mobility Cup SATURNE® II - Size 42	1-0111342
Dual Mobility Cup SATURNE® II - Size 44	1-0111344
Dual Mobility Cup SATURNE® II - Size 46	1-0111346
Dual Mobility Cup SATURNE® II - Size 48	1-0111348
Dual Mobility Cup SATURNE® II - Size 50	1-0111350
Dual Mobility Cup SATURNE® II - Size 52	1-0111352
Dual Mobility Cup SATURNE® II - Size 54	1-0111354
Dual Mobility Cup SATURNE® II - Size 56	1-0111356
Dual Mobility Cup SATURNE® II - Size 58	1-0111358
Dual Mobility Cup SATURNE® II - Size 60	1-0111360
Dual Mobility Cup SATURNE® II - Size 62	1-0111362
Dual Mobility Cup SATURNE® II - Size 64	1-0111364

DUAL MOBILITY LINERS SATURNE®

	Femoral heads Ø22,2mm	Femoral heads Ø28mm
Liner SATURNE® - Size 42	1-0103642	-
Liner SATURNE® - Size 44	1-0103644	-
Liner SATURNE® - Size 46	1-0103646	-
Liner SATURNE® - Size 48	1-0103648	1-0103748
Liner SATURNE® - Size 50	1-0103650	1-0103750
Liner SATURNE® - Size 52	1-0103652	1-0103752
Liner SATURNE® - Size 54	1-0103654	1-0103754
Liner SATURNE® - Size 56	1-0103656	1-0103756
Liner SATURNE® - Size 58	1-0103658	1-0103758
Liner SATURNE® - Size 60	1-0103660	1-0103760
Liner SATURNE® - Size 62	1-0103662	1-0103762
Liner SATURNE® - Size 64	1-0103664	1-0103764

Creation date : December 2016—Legal manufacturer : AMPLITUDE, 11 cours Jacques Offenbach, 26000 Valence, FRANCE—Destination : replacement of hip joint in association with a femoral component

Customer Service - France

Porte du Grand Lyon,
01700 Neyron, France
Tel. : +33 (0)4 37 85 19 19
Fax : +33 (0)4 37 85 19 18
E-mail : amplitude@amplitude-ortho.com

Customer Service-Export :

11, cours Jacques Offenbach. Zone Mozart 2,
26000 Valence, France
Tel. : +33 (0)4 75 41 87 41
Fax : +33 (0)4 75 41 87 42
E-mail : amplitude@amplitude-ortho.com

Internet : www.amplitude-ortho.com

 AMPLITUDE®