



TECHNICAL SHEET

 **horizon II**




Hemispherical acetabular cup with **flat pole**

Equatorial press fit

3 to 5 holes to ensure **primary fixation** with screws


Cementless Cup



9 sizes available

Material: Anodized Titanium Alloy (Ti6Al4V)

Important Features



Fixation Screws



Fixation screws $\varnothing$ 6.5mm

Use of 1 to 5 screws

Material: Anodized Titanium alloy (Ti6Al4V)

HORIZON II with holes

HORIZON II Acetabular Cup with Holes - TI/HA Cementless - Size 48	1-01102 48
HORIZON II Acetabular Cup with Holes - TI/HA Cementless - Size 50	1-01102 50
HORIZON II Acetabular Cup with Holes - TI/HA Cementless - Size 52	1-01102 52
HORIZON II Acetabular Cup with Holes - TI/HA Cementless - Size 54	1-01102 54
HORIZON II Acetabular Cup with Holes - TI/HA Cementless - Size 56	1-01102 56
HORIZON II Acetabular Cup with Holes - TI/HA Cementless - Size 58	1-01102 58
HORIZON II Acetabular Cup with Holes - TI/HA Cementless - Size 60	1-01102 60
HORIZON II Acetabular Cup with Holes - TI/HA Cementless - Size 62	1-01102 62
HORIZON II Acetabular Cup with Holes - TI/HA Cementless - Size 64*	1-01102 64*

BIOLOX [®] Delta ceramic liner	Ø 28	Ø 32	Ø 36
BIOLOX [®] Delta Ceramic Liner for HORIZON II Acetabular Cup - Size 48	1-0108 636	-	-
BIOLOX [®] Delta Ceramic Liner for HORIZON II Acetabular Cup - Size 50 and 52	-	1-01100 39	-
BIOLOX [®] Delta Ceramic Liner for HORIZON II Acetabular Cup - Size 54, 56 and 58	-	1-01100 44	1-01101 44
BIOLOX [®] Delta Ceramic Liner for HORIZON II Acetabular Cup - Size 60 and 62	-	-	1-01101 48
BIOLOX [®] Delta Ceramic Liner for HORIZON II Acetabular Cup - Size 64*	-	-	1-01101 52*

Fixation screws

Fixation screw Ø 6.5 mm - Lg 16 mm	1-01016 16
Fixation screw Ø 6.5 mm - Lg 20 mm	1-01016 20
Fixation screw Ø 6.5 mm - Lg 25 mm	1-01016 25
Fixation screw Ø 6.5 mm - Lg 30 mm	1-01016 30
Fixation screw Ø 6.5 mm - Lg 35 mm	1-01016 35
Fixation screw Ø 6.5 mm - Lg 40 mm	1-01016 40
Fixation screw Ø 6.5 mm - Lg 45 mm	1-01016 45

**optional, available on request only*

Service Clients – France :

Porte du Grand Lyon,
01700 Neyron – France
Tél. : **+33 (0)4 37 85 19 19**
Fax : **+33 (0)4 37 85 19 18**
E-mail : amplitude@amplitude-ortho.com

Customer Service – Export :

11 cours Jacques Offenbach,
ZA Mozart 2,
26000 Valence – France
Tél. : **+33 (0)4 75 41 87 41**
Fax : **+33 (0)4 75 41 87 42**

BIOLOX[®] is a registered trademark of the CeramTec Group

www.amplitude-ortho.com