

SCORE® AS (Allergy Solution)



- SCORE® AS Total knee prosthesis coated with Titanium Nitride (TiN) ceramic to minimise allergic reactions
- Coating acts as a barrier to reduce the release of allergenic ions
- Coating thickness: 4 µm
- Cemented components:
 - Femoral component
 - Tibial base plate
 - Tibial extension stems
 - Tibial offset connectors
 - Tibial half blocks
- Tibial offset connectors allows the lateralization of the stem in the space between the diaphyseal axis and the centre of the knee
- Highly polished articular surfaces:
 - Femoral component / tibial insert
 - Tibial Insert / tibial base plate
- Available only upon request

SCORE® Allergy Solution femoral component PCL sacrificing cemented Size 1 to Size 7 RIGHT	1-0206101 to 1-0206107
SCORE® Allergy Solution femoral component PCL sacrificing cemented Size 1 to Size 7 LEFT	1-0206111 to 1-0206117
SCORE® Allergy Solution tibial base plate for mobile bearing insert cemented Size 1 to Size 7	1-0206201 to 1-0206207
SCORE® mobile bearing insert - Size 1 to 7 thickness 10 mm	1-0200611 to 1-0200617
SCORE® mobile bearing insert - Size 1 to 7 thickness 12 mm	1-0200621 to 1-0200627
SCORE® mobile bearing insert - Size 1 to 7 thickness 14 mm	1-0200631 to 1-0200637
SCORE® mobile bearing insert - Size 1 to 7 thickness 16 mm	1-0200641 to 1-0200641
SCORE® mobile bearing insert - Size 1 to 7 thickness 20 mm	1-0200651 to 1-0200651

Available upon request

Allergy Solution extension stem Ø 10 length 75 mm	1-0207310
Allergy Solution extension stem Ø 10 length 100 mm	1-0207311
Allergy Solution extension stem Ø 10 length 150 mm	1-0207312
Allergy Solution extension stem Ø 12 length 75 mm	1-0207320
Allergy Solution extension stem Ø 12 length 100 mm	1-0207321
Allergy Solution extension stem Ø 12 length 150 mm	1-0207322
Allergy Solution extension stem Ø 12 length 200 mm	1-0207323
Allergy Solution extension stem Ø 14 length 75 mm	1-0207330
Allergy Solution extension stem Ø 14 length 100 mm	1-0207331
Allergy Solution extension stem Ø 14 length 150 mm	1-0207332
Allergy Solution extension stem Ø 14 length 200 mm	1-0207333
Allergy Solution extension stem Ø 16 length 100 mm	1-0207341
Allergy Solution extension stem Ø 16 length 150 mm	1-0207342
Allergy Solution extension stem Ø 16 length 200 mm	1-0207343
Allergy Solution offset connector 2 mm	1-0207402
Allergy Solution offset connector 4 mm	1-0207404
Allergy Solution offset connector 6 mm	1-0207406
Allergy Solution tibial half block Size 1/2 thickness 5 mm	1-0207511
Allergy Solution tibial half block Size 3/4/5 thickness 5 mm	1-0207513
Allergy Solution tibial half block Size 6/7 thickness 5 mm	1-0207516
Allergy Solution tibial half block Size 1/2 thickness 10 mm MR/LL	1-0207521
Allergy Solution tibial half block Size 1/2 thickness 15 mm MR/LL	1-0207531
Allergy Solution tibial half block Size 3/4/5 thickness 10 mm MR/LL	1-0207523
Allergy Solution tibial half block Size 3/4/5 thickness 15 mm MR/LL	1-0207533
Allergy Solution tibial half block Size 6/7 thickness 10 mm MR/LL	1-0207526
Allergy Solution tibial half block Size 6/7 thickness 15 mm MR/LL	1-0207536
Allergy Solution tibial half block Size 1/2 thickness 10 mm LR/ML	1-0207541
Allergy Solution tibial half block Size 1/2 thickness 15 mm LR/ML	1-0207551
Allergy Solution tibial half block Size 3/4/5 thickness 10 mm LR/ML	1-0207543
Allergy Solution tibial half block Size 3/4/5 thickness 15 mm LR/ML	1-0207553
Allergy Solution tibial half block Size 6/7 thickness 10 mm LR/ML	1-0207546
Allergy Solution tibial half block Size 6/7 thickness 15 mm LR/ML	1-0207556