



SPHERIC® RL : features



- Bipolar cup for hemiarthroplasty of the hip
- Outer surface: highly polished
- Pre-assembled constrained liner for $\varnothing 28$ mm femoral heads
- 9 sizes (odd numbers): 43 to 59 mm
- Material: Stainless steel (316L) and ultra-high molecular weight polyethylene (UHMWPE)

Reference : FT.GB028/1.0

SPHERIC® RL : implants

CE 0499

SPHERIC® RL Constrained Bipolar Cup — Size 43/28	1-0191543
SPHERIC® RL Constrained Bipolar Cup — Size 45/28	1-0191545
SPHERIC® RL Constrained Bipolar Cup — Size 47/28	1-0191547
SPHERIC® RL Constrained Bipolar Cup — Size 49/28	1-0191549
SPHERIC® RL Constrained Bipolar Cup — Size 51/28	1-0191551
SPHERIC® RL Constrained Bipolar Cup — Size 53/28	1-0191553
SPHERIC® RL Constrained Bipolar Cup — Size 55/28	1-0191555
SPHERIC® RL Constrained Bipolar Cup — Size 57/28	1-0191557
SPHERIC® RL Constrained Bipolar Cup — Size 59/28	1-0191559

Manufacturer: EVOLUTIS, av. de la Libération, 42720 Briennon, France