




AMPLITUDE[®]

GENERIC[®]
femoral stem
cemented
smooth

GENERIC® : cemented femoral stem




- Neck-shaft angle 137°
- Highly polished surface
- Material : M30NW

GENERIC® : main elements



- Mirror polished thin neck



- Small-sized taper (10/12)5°42'30"



- Threaded hole for extraction

GENERIC® femoral stem - cimented - smooth - size 1	1-0100201
GENERIC® femoral stem - cimented - smooth - size 2	1-0100202
GENERIC® femoral stem - cimented - smooth - size 3	1-0100203
GENERIC® femoral stem - cimented - smooth - size 4	1-0100204
GENERIC® femoral stem - cimented - smooth - size 5	1-0100205
GENERIC® femoral stem - cimented - smooth - size 6	1-0100206
GENERIC® femoral stem - cimented - smooth - size 7	1-0100207

M30NW metallic femoral head - Ø 22,2 short neck	1-0101101
M30NW metallic femoral head - Ø 22,2 medium neck	1-0101102
M30NW metallic femoral head - Ø 22,2 long neck	1-0101103
M30NW metallic femoral head - Ø 28 short neck	1-0101111
M30NW metallic femoral head - Ø 28 medium neck	1-0101112
M30NW metallic femoral head - Ø 28 long neck	1-0101113
M30NW metallic femoral head - Ø 28 extra long neck	1-0101114
M30NW metallic femoral head - Ø 32 short neck	1-0101121
M30NW metallic femoral head - Ø 32 medium neck	1-0101122
M30NW metallic femoral head - Ø 32 long neck	1-0101123

Delta ceramic femoral head - Ø 28 Col short neck	1-0191101
Delta ceramic femoral head - Ø 28 Col medium neck	1-0191102
Delta ceramic femoral head - Ø 28 Col long neck	1-0191103
Delta ceramic femoral head - Ø 32 Col short neck	1-0191111
Delta ceramic femoral head - Ø 32 Col medium neck	1-0191122
Delta ceramic femoral head - Ø 32 Col long neck	1-0191113
Delta ceramic femoral head - Ø 36 Col short neck	1-0191121
Delta ceramic femoral head - Ø 36 Col medium neck	1-0191122
Delta ceramic femoral head - Ø 36 Col long neck	1-0191123