




AMPLITUDE ©

**SATURNE®**  
Reconstruction **SATURNE®**  
Double mobility cups  
Cementless version



## SATURNE® cup: main elements

---



- Anti-dislocation anatomical shape
- Flattened pole
- Equatorial press-fit
- Double coating: Titanium spray + HAP (80 µm)



- Mirror-polished inner surface
- Rounded edges
- Material : M30NW



- Liner with a chamfer:
  - to optimise the range of motion
  - to limit wear to a minimum
- Retentive liner for  $\varnothing$  22.2 and 28 mm heads



## Reconstruction SATURNE® cup: main elements

---



- Anti-dislocation anatomical shape
- Flattened pole
- Equatorial press-fit
- Double coating: Titanium spray + HAP (80 µm)
- Mirror-polished inner surface
- Rounded edges
- Material : M25NW



- On the upper acetabular rim, 2 flanges:
  - to be cut and bent
  - located at the equator of the shell for easier positioning
- With two holes for fixation with the use of  $\varnothing$  4,5 mm cortical screws
- In the inferior acetabular rim, a positioning hook located in the obturator hole
- Liner with a chamfer:
  - to optimise the range of motion
  - to limit wear to a minimum
- Retentive liner for  $\varnothing$  22.2 and 28 mm heads

## SATURNE® and reconstruction SATURNE® cups: instrumentation C€

---



- Navigated instrumentation



## SATURNE® cup: implants

CE 0086

SATURNE® Press- fit Double mobility cup size 44	1-0103444	Liner for SATURNE® Double mobility cup size 44/22	1-0103644
SATURNE® Press- fit Double mobility cup size 46	1-0103446	Liner for SATURNE® Double mobility cup size 46/22	1-0103646
SATURNE® Press- fit Double mobility cup size 48	1-0103448	Liner for SATURNE® Double mobility cup size 48/22	1-0103648
SATURNE® Press- fit Double mobility cup size 50	1-0103450	Liner for SATURNE® Double mobility cup size 50/22	1-0103650
SATURNE® Press- fit Double mobility cup size 52	1-0103452	Liner for SATURNE® Double mobility cup size 52/22	1-0103652
SATURNE® Press- fit Double mobility cup size 54	1-0103454	Liner for SATURNE® Double mobility cup size 54/22	1-0103654
SATURNE® Press- fit Double mobility cup size 56	1-0103456	Liner for SATURNE® Double mobility cup size 56/22	1-0103656
SATURNE® Press- fit Double mobility cup size 58	1-0103458	Liner for SATURNE® Double mobility cup size 58/22	1-0103658
SATURNE® Press- fit Double mobility cup size 60	1-0103460	Liner for SATURNE® Double mobility cup size 60/22	1-0103660
SATURNE® Press- fit Double mobility cup size 62	1-0103462	Liner for SATURNE® Double mobility cup size 62/22	1-0103662
SATURNE® Press- fit Double mobility cup size 64	1-0103464	Liner for SATURNE® Double mobility cup size 64/22	1-0103664
		Liner for SATURNE® Double mobility cup size 48/28	1-0103748
		Liner for SATURNE® Double mobility cup size 50/28	1-0103750
		Liner for SATURNE® Double mobility cup size 52/28	1-0103752
		Liner for SATURNE® Double mobility cup size 54/28	1-0103754
		Liner for SATURNE® Double mobility cup size 56/28	1-0103756
		Liner for SATURNE® Double mobility cup size 58/28	1-0103758
		Liner for SATURNE® Double mobility cup size 60/28	1-0103760
		Liner for SATURNE® Double mobility cup size 62/28	1-0103762
		Liner for SATURNE® Double mobility cup size 64/28	1-0103764

## Reconstruction SATURNE® cup: implants

CE 0086

Reconstruction SATURNE® Double mobility cup size 50	1-0103350	Liner for SATURNE® Double mobility cup size 50/22	1-0103650
Reconstruction SATURNE® Double mobility cup size 52	1-0103352	Liner for SATURNE® Double mobility cup size 52/22	1-0103652
Reconstruction SATURNE® Double mobility cup size 54	1-0103354	Liner for SATURNE® Double mobility cup size 54/22	1-0103654
Reconstruction SATURNE® Double mobility cup size 56	1-0103356	Liner for SATURNE® Double mobility cup size 56/22	1-0103656
Reconstruction SATURNE® Double mobility cup size 58	1-0103358	Liner for SATURNE® Double mobility cup size 58/22	1-0103658
Reconstruction SATURNE® Double mobility cup size 60	1-0103360	Liner for SATURNE® Double mobility cup size 60/22	1-0103660
Reconstruction SATURNE® Double mobility cup size 62	1-0103362	Liner for SATURNE® Double mobility cup size 62/22	1-0103662
Reconstruction SATURNE® Double mobility cup size 64	1-0103364	Liner for SATURNE® Double mobility cup size 64/22	1-0103664
Reconstruction SATURNE® Double mobility cup size 66	1-0103366	Liner for SATURNE® Double mobility cup size 66/22	1-0103666
		Liner for SATURNE® Double mobility cup size 50/28	1-0103750
		Liner for SATURNE® Double mobility cup size 52/28	1-0103752
		Liner for SATURNE® Double mobility cup size 54/28	1-0103754
		Liner for SATURNE® Double mobility cup size 56/28	1-0103756
		Liner for SATURNE® Double mobility cup size 58/28	1-0103758
		Liner for SATURNE® Double mobility cup size 60/28	1-0103760
		Liner for SATURNE® Double mobility cup size 62/28	1-0103762
		Liner for SATURNE® Double mobility cup size 64/28	1-0103764
		Liner for SATURNE® Double mobility cup size 66/28	1-0103766

Distributed by :



**AMPLITUDE** ©

Customer Service France : Porte du Grand Lyon, 01700 Neyron, France. Tél. : +33 (0)4 37 85 19 19. Fax : +33 (0)4 37 85 19 18

Customer Service Export : 11, cours Jacques Offenbach, Zone Mozart 2, 26000 Valence, France. Tél. : +33 (0)4 75 41 87 41. Fax : +33 (0)4 75 41 87 42

E-mail : [amplitude@amplitude-ortho.com](mailto:amplitude@amplitude-ortho.com)

Internet : [www.amplitude-ortho.com](http://www.amplitude-ortho.com)