




AMPLITUDE ©

SATURNE®
Cemented tripolar cup

SATURNE® acetabular component and acetabular ring



- Anatomical anti-dislocation profile.
- Flattened pole.
- Rounded up edge of the cup.
- Grooves for a homogeneous distribution of the cement.
- M30NW cup : mirror-polish inner surface and sand-blast exterior surface.
- No hole.



- Insert for \varnothing 22,2 mm or \varnothing 28 mm femoral head (metallic or alumina).
- Peripheral bevelled protecting the insert's retention area.
- UHMWPE insert.



- Anatomic acetabular rings : left and right
- Anterior branch shorter than the posterior one.
- M30NW hemispherical cross 2 mm thick, with :
 - in the lower part a positioning foramen hook.
 - in the upper part a plate with 4 holes for fixation with 4,5 mm cortical screw.
- Sterile packaging.


SATURNE® acetabular component



- The cemented SATURNE® cup combined with the acetabular ring enables the dual mobility for bone reconstruction with requires a cement fixation.

SATURNE® acetabular component: instrumentation

CE



- Navigable instrumentation identical for cemented and cementless SATURNE® cups and for the reconstruction SATURNE® cup.
- Minimal obstruction for the reduced access ways.

Acetabular ring: instrumentation

CE



- Instrumentation allowing the confirmation of the size and the orientation of the implant.

SATURNE® acetabular component: implants

CE 0086

SATURNE® tripolar cup cemented - size 44	1-0106244
SATURNE® tripolar cup cemented - size 46	1-0106246
SATURNE® tripolar cup cemented - size 48	1-0106248
SATURNE® tripolar cup cemented - size 50	1-0106250
SATURNE® tripolar cup cemented - size 52	1-0106252
SATURNE® tripolar cup cemented - size 54	1-0106254
SATURNE® tripolar cup cemented - size 56	1-0106256
SATURNE® tripolar cup cemented - size 58	1-0106258
SATURNE® tripolar cup cemented - size 60	1-0106260

SATURNE® insert for tripolar cup - size 44/22	1-0103644
SATURNE® insert for tripolar cup - size 46/22	1-0103646
SATURNE® insert for tripolar cup - size 48/22	1-0103648
SATURNE® insert for tripolar cup - size 50/22	1-0103650
SATURNE® insert for tripolar cup - size 52/22	1-0103652
SATURNE® insert for tripolar cup - size 54/22	1-0103654
SATURNE® insert for tripolar cup - size 56/22	1-0103656
SATURNE® insert for tripolar cup - size 58/22	1-0103658
SATURNE® insert for tripolar cup - size 60/22	1-0103660

SATURNE® insert for tripolar cup - size 48/28	1-0103748
SATURNE® insert for tripolar cup - size 50/28	1-0103750
SATURNE® insert for tripolar cup - size 52/28	1-0103752
SATURNE® insert for tripolar cup - size 54/28	1-0103754
SATURNE® insert for tripolar cup - size 56/28	1-0103756
SATURNE® insert for tripolar cup - size 58/28	1-0103758
SATURNE® insert for tripolar cup - size 60/28	1-0103760

Acetabular ring: implants

CE 0086

Right acetabular ring for cup size 37/40	1-0103140
Right acetabular ring for cup size 40/44	1-0103144
Right acetabular ring for cup size 43/48	1-0103148
Right acetabular ring for cup size 46/52	1-0103152
Right acetabular ring for cup size 49/56	1-0103156
Right acetabular ring for cup size 53/60	1-0103160

Left acetabular ring for cup size 37/40	1-0103240
Left acetabular ring for cup size 40/44	1-0103244
Left acetabular ring for cup size 43/48	1-0103248
Left acetabular ring for cup size 46/52	1-0103252
Left acetabular ring for cup size 49/56	1-0103256
Left acetabular ring for cup size 53/60	1-0103260

Distributed by:

