





AMPLITUDE[®]


LOGIC
femoral stem
cementless
partial HA coated

LOGIC® femoral stem cementless partial HA coated



- Neck-shaft angle 137°
- Highly polished distal stem
- Blind proximal grooves
- 80 µm HAP coating in the metaphyseal region
- 7 sizes : homothetic increase in sizes
- Material : anodized TA6V

LOGIC® : Main elements



- Mirror polished thin neck



- Small-sized taper (10/12) 5-42°



- Threaded hole for extraction

LOGIC® femoral stem partial HA coated - cementless - size 1	1-0100101
LOGIC® femoral stem partial HA coated - cementless - size 2	1-0106502
LOGIC® femoral stem partial HA coated - cementless - size 3	1-0106503
LOGIC® femoral stem partial HA coated - cementless - size 4	1-0106504
LOGIC® femoral stem partial HA coated - cementless - size 5	1-0106505
LOGIC® femoral stem partial HA coated - cementless - size 6	1-0106506
LOGIC® femoral stem partial HA coated - cementless - size 7	1-0106507

M30NW metallic femoral head - Ø 22,2 short neck	1-0101101
M30NW metallic femoral head - Ø 22,2 medium neck	1-0101102
M30NW metallic femoral head - Ø 22,2 long neck	1-0101103
M30NW metallic femoral head - Ø 28 short neck	1-0101111
M30NW metallic femoral head - Ø 28 medium neck	1-0101112
M30NW metallic femoral head - Ø 28 long neck	1-0101113
M30NW metallic femoral head - Ø 28 extra long neck	1-0101114
M30NW metallic femoral head - Ø 32 short neck	1-0101121
M30NW metallic femoral head - Ø 32 medium neck	1-0101122
M30NW metallic femoral head - Ø 32 long neck	1-0101123

Delta ceramic femoral head - Ø 28 Col short neck	1-0191101
Delta ceramic femoral head - Ø 28 Col medium neck	1-0191102
Delta ceramic femoral head - Ø 28 Col long neck	1-0191103
Delta ceramic femoral head - Ø 32 Col short neck	1-0191111
Delta ceramic femoral head - Ø 32 Col medium neck	1-0191122
Delta ceramic femoral head - Ø 32 Col long neck	1-0191113
Delta ceramic femoral head - Ø 36 Col short neck	1-0191121
Delta ceramic femoral head - Ø 36 Col medium neck	1-0191122
Delta ceramic femoral head - Ø 36 Col long neck	1-0191123