

Technical sheet

# SPHERIC<sup>®</sup>

---

Bipolar cup



  
AMPLITUDE<sup>®</sup>

# References

- Bipolar cup for hip hemiarthroplasties
- 2 articulation levels:
  - First articulation of the femoral head (22,2 mm or 28mm) within the preassembled liner
  - Second articulation of the external mirror-polished surface, freely in the acetabulum: the articular diameter is then as large as possible for maximum joint stability
- Retentivity of the head is ensured by a retention ring, assembled during surgery, at the threshold of the head cavity
- Materials:
  - Metallic cup: Stainless steel (M30NW)
  - Liner: Ultra-high molecular weight polyethylene (UHMWPE)



CE<sub>2797</sub>

Reference : FT.CB.024/1.1

## SPHERIC®

	Ø22,2 mm	Ø28 mm
SPHERIC® bipolar cup - Size <b>40</b>	1-01 <b>02540</b>	-
SPHERIC® bipolar cup - Size <b>42</b>	1-01 <b>02542</b>	-
SPHERIC® bipolar cup - Size <b>44</b>	1-01 <b>02544</b>	-
SPHERIC® bipolar cup - Size <b>46</b>	1-01 <b>02546</b>	-
SPHERIC® bipolar cup - Size <b>48</b>	1-01 <b>02548</b>	1-01 <b>02648</b>
SPHERIC® bipolar cup - Size <b>50</b>	1-01 <b>02550</b>	1-01 <b>02650</b>
SPHERIC® bipolar cup - Size <b>52</b>	1-01 <b>02552</b>	1-01 <b>02652</b>
SPHERIC® bipolar cup - Size <b>54</b>	1-01 <b>02554</b>	1-01 <b>02654</b>
SPHERIC® bipolar cup - Size <b>56</b>	1-01 <b>02556</b>	1-01 <b>02656</b>
SPHERIC® bipolar cup - Size <b>58</b>	-	1-01 <b>02658</b>

**Service Clients-France :**  
 Porte du Grand Lyon  
 01700 Neyron, France  
 Tel. : +33 (0)4 37 85 19 19  
 Fax : +33 (0)4 37 85 19 18  
 E-mail : [amplitude@amplitude-ortho.com](mailto:amplitude@amplitude-ortho.com)

**Customer Service-Export :**  
 11, cours Jacques Offenbach  
 Zone Mozart 2  
 26000 Valence, France  
 Tel. : +33 (0)4 75 41 87 41  
 Fax : +33 (0)4 75 41 87 42  
 E-mail : [amplitude@amplitude-ortho.com](mailto:amplitude@amplitude-ortho.com)

Internet : [www.amplitude-ortho.com](http://www.amplitude-ortho.com)