

ECLIPSE[®] Profil

Resorbable polylactic interference screw




Distributed by

 **AMPLITUDE[®]**

Features

- Canulated Polylactic interference screw
- screw with
 - Progressive insert
 - Automatic head
 - Pentalobe cannulated



Material

PLDLLA

- PLDLLA, polylactic acide L & D
- “ L-Lactide provides greater mechanical resistance. DL-Lactide is amorphous (the crystallinity has been removed) which increases the resorption rate. Biocompatibility improvement (amorphous D,L-lactide units) compared to a pure PLLA


Range :

- Diameters 7/8/9 mm Length 20/25 and 30mm
- Diameters 10mm Length 25 et 30mm
- Diameters 11 et 12mm Length 35mm

Diamètres	Ø7mm	Ø8mm	Ø9mm	Ø10mm	Ø11mm	Ø12mm
Longueurs						
20 mm						
Références	1-0402771	1-0402781	1-0402791			
25mm						
Références	1-0402772	1-0402782	1-0402792	1-0402702		
30mm						
Références	1-0402773	1-0402783	1-0402793	1-0402703		
35mm						
Références					1-0402714	1-0402724

Screwing instrumentation : the choice

- Ratcheting trinkle driver handle : Amplitude
 - Shank screwdriver
 - Shank starter



- Screw Guidewire 1.1 mm diameter
 - Rigid : Amplitude
 - Nitinol : Biomatlante

- Monobloc: Biomatlante
 - screwdriver
 - starter



References

ECLIPSE® Profil " Resorbable polylactic interference screw

Manufactured by Biomatlante 

Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 7 Lg 20 mm	1-0402771
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 7 Lg 25 mm	1-0402772
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 7 Lg 30 mm	1-0402773
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 8 Lg 20 mm	1-0402781
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 8 Lg 25 mm	1-0402782
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 8 Lg 30 mm	1-0402783
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 9 Lg 20 mm	1-0402791
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 9 Lg 25 mm	1-0402792
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 9 Lg 30 mm	1-0402793
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 10 Lg 25 mm	1-0402702
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 10 Lg 30 mm	1-0402703
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 11 Lg 35 mm	1-0402714
Resorbable Polylactique Interference screw ECLIPSE® Profil Ø 12 Lg 35 mm	1-0402724

Instrumentation Set

Manufactured by AMPLITUDE 

Modular Instrumentation Set In-Out	2-0299940
Modular Instrumentation Set Out-In	2-0299957

Creation date : March 2019 – Composite interference screw BCP Legal manufacturer : Biomatlante • 5 rue Édouard Belin • Z.A. Les IV Nations • 44360 • Vigneux-de-Bretagne • France – Destination : medical device for Knee ligament surgery
–instrumentation set legal manufacturers : Amplitude , 11 cours Jacques Offenbach, 26000 Valence, France and Biomatlante • 5 rue Édouard Belin • Z.A. Les IV Nations • 44360 • Vigneux-de-Bretagne • France

Customer Service-France :

Porte du Grand Lyon,
01700 Neyron, France
Tel. : +33 (0)4 37 85 19 19
Fax : +33 (0)4 37 85 19 18

E-mail : amplitude@amplitude-ortho.com

Customer Service-Export :

11, cours Jacques Offenbach. Zone Mozart 2,
26000 Valence, France
Tel. : +33 (0)4 75 41 87 41
Fax : +33 (0)4 75 41 87 42

Internet : www.amplitude-ortho.com