



AMPLITUDE<sup>©</sup>


**GENERIC<sup>®</sup> Revision**  
cemented  
**INTEGRALE<sup>®</sup> Revision**  
cementless

Femoral revision stems




## Main elements


---



- Neck-shaft angle  $135^{\circ}$
- Perfect visualization of the impact limit
- Stainless steel (M30NW) or ceramic femoral head



- 7 stem sizes with homothetically increasing dimensions
- Small sized proximal taper: 10/12 taper ( $5^{\circ}42'30''$ )
- Narrow mirror-like polished neck design



- Threaded hole for easier extraction
- Ovoid metaphyseal section



## GENERIC® revision and INTEGRALE® revision stems

---



### GENERIC® Revision:

- Highly polished surface
- Material: M30NW

### INTEGRALE® Revision:

- Specific grooves to increase the bone contact surface
- Material: anodized TA6V
- HAP coating (80 µm)



## GENERIC® revision and INTEGRALE® revision stems: implants CE 0086

Cemented Revision GENERIC® modular femoral stem Size 1	1-0100401
Cemented Revision GENERIC® modular femoral stem Size 2	1-0100402
Cemented Revision GENERIC® modular femoral stem Size 3	1-0100403
Cemented Revision GENERIC® modular femoral stem Size 4	1-0100404
Cemented Revision GENERIC® modular femoral stem Size 5	1-0100405
Cemented Revision GENERIC® modular femoral stem Size 6	1-0100406
Cemented Revision GENERIC® modular femoral stem Size 7	1-0100407

Cementless HAP Revision INTEGRALE® modular femoral stem Size 1	1-0100301
Cementless HAP Revision INTEGRALE® modular femoral stem Size 2	1-0100302
Cementless HAP Revision INTEGRALE® modular femoral stem Size 3	1-0100303
Cementless HAP Revision INTEGRALE® modular femoral stem Size 4	1-0100304
Cementless HAP Revision INTEGRALE® modular femoral stem Size 5	1-0100305
Cementless HAP Revision INTEGRALE® modular femoral stem Size 6	1-0100306
Cementless HAP Revision INTEGRALE® modular femoral stem Size 7	1-0100307

Metallic femoral head M30NW Ø 22.2 Short Neck	1-0101101
Metallic femoral head M30NW Ø 22.2 Medium Neck	1-0101102
Metallic femoral head M30NW Ø 22.2 Long Neck	1-0101103
Metallic femoral head M30NW Ø 28 Short Neck	1-0101111
Metallic femoral head M30NW Ø 28 Medium Neck	1-0101112
Metallic femoral head M30NW Ø 28 Long Neck	1-0101113
Metallic femoral head M30NW Ø 28 Extra Long Neck	1-0101114
Metallic femoral head M30NW Ø 32 Short Neck	1-0101121
Metallic femoral head M30NW Ø 32 Medium Neck	1-0101122
Metallic femoral head M30NW Ø 32 Long Neck	1-0101123

Ceramic femoral head Ø 28 Short Neck	1-0102711
Ceramic femoral head Ø 28 Medium Neck	1-0102712
Ceramic femoral head Ø 28 Long Neck	1-0102713
Ceramic femoral head Ø 32 Short Neck	1-0102721
Ceramic femoral head Ø 32 Medium Neck	1-0102722
Ceramic femoral head Ø 32 Long Neck	1-0102723
Ceramic femoral head Ø 36 Short Neck	1-0102731
Ceramic femoral head Ø 36 Medium Neck	1-0102732
Ceramic femoral head Ø 36 Long Neck	1-0102733

Distributed by:



**AMPLITUDE** 

Customer Service France: Porte du Grand Lyon, 01700 Neyron, France. Tel.: +33 (0)4 37 85 19 19. Fax: +33 (0)4 37 85 19 18

Customer Service Export: 11, cours Jacques Offenbach, Zone Mozart 2, 26000 Valence, France. Tel.: +33 (0)4 74 41 87 41. Fax: +33 (0)4 75 41 87 42

E-mail: [amplitude@amplitude-ortho.com](mailto:amplitude@amplitude-ortho.com)

Internet: [www.amplitude-ortho.com](http://www.amplitude-ortho.com)